




Packaging instructions TWINWARM BB





Status: Juni 2024



Based on work instruction BB5400de_rev.03

DO507 Rev. 02

Goal	Packaging of the TWINWARM BB device for transportation with Instapack foam bags
Work equipment	Packaging tape Packaging tape hand dispenser
a) Components	TWINWARM BB shipping box (dimensions: 500 x 500 x 550 mm) (item no. V000034) Instapack foam bag (3 pieces) (art. no. V000038)
b) Tools	1. optional: side cutter
General safety instructions	 Safety
	 Note
	 Quality / Execution

Work process

No.	Process step	Description	Schematic / Photo
1	Prepare TWINWARM BB	<p>+ Ensure that the mains plug has been disconnected.</p>	
		Fill in and print out the repair order.	
2	Foam bag placement	The Instapack foam bags are placed in the box around the appliance as follows.	
3	Foam bag placement	1) Place the bottom cushion in the box.	
		2) Place the appliance in the box. (The rollers fit into the recesses).	
		3) Then place the top cushion on the appliance.	

4	Additional accessories	Attach the repair order	
5	Seal the box	Seal the box with suitable cardboard tape.	
6	Examination	<p>◆ Visually check the condition of the appliance packaging. Take a photo as proof in the event of any transport damage.</p>	