

Warming blankets and mats

temperature control from above and below with pressure-relief

INFANT CARDIO

Type	– 909
Size* (LxW)	– Thorax hoses 48 x 19 cm
	– Lower section 86 x 73 cm
Weight	– about 590 g



- The infant cardio mat of the **MoeckWarmingSystem**[®] is fitted with three air inlets
- Designed for the special requirements of cardiac catheter tables, this mat contains lateral rolls of fabric to prevent liquid spill. It is also fitted with a reclining area with perforations for liquid drainage
- The integrated thorax hoses can be draped around the patient at any time. They offer even temperature control from above or from the side
- The flaps on the side of the mat can either be used to fix it onto the OR table or simply be turned down.

Tip: If the mat is turned 180°, the thorax hoses can also be used for temperature control of the lower extremities.



* since the textile products shrink during the first specialized washing cycles, the specified dimensions of the usable surface (excluding air inlets) apply after several cycles and with a tolerance of ± 2 cm.