






Description	Warming mat for temperature control from below
Type No.	906
Art. No.	11WD906M
Dimensions¹	204×64 cm (L×W)
Weight	approx. 1,960 g
Materials	Green Textile and padding: 100% polyester / White Textile: 100% polyester with polyurethane coating
Hose connections	● ● ●
Features	► Placement under the patient ► Full, unrestricted access to the patient ► The integrated layer reduces contact pressure and minimises fluid accumulation on the surface

Examples of use²



Application notes³

Pre-OP	Intra-OP	Post-OP	Supine Position	Prone Position	Lateral Position	Lithotomy Position	Beachchair Position	Transfer	Paediatrics
General surgery			✓		Neurosurgery			✓	
Vascular surgery			✓		Orthopaedics and Trauma surgery			✓	
Obstetrics and Gynaecology			✓		Paediatrics				
Cardiac surgery			✓		Plastic surgery			✓	
Head surgery			✓		Abdominal surgery			✓	
Maxillofacial surgery			✓		Thoracic surgery			✓	



¹ Fabric products shrink by 3–5% during the first reprocessing cycle. A tolerance of ±2 cm also applies to the stated dimensions.

² The examples of use show all possible connections available — one to max. two hoses are sufficient for proper operation.

³ The application notes are based on the practical experience of our service partners and are therefore not binding.